

Aerosols into Suspension Collectors, a New Approach for and Efficient Collection of Airborne Particles for Toxicological Studies

Arantzazu Eiguren Fernandez¹, Nicholas Buglak², Steven Spielman¹,
Gregory Lewis¹, Ning Li² and Susanne Hering¹

¹Aerosol Dynamics Inc., Berkeley, CA

²Department of Pathology and Diagnostic Investigation, Michigan State University, MI

BACKGROUND


- Accurate information on the toxicological responses to ambient aerosols is needed for human exposure models and risk assessment studies
- In vitro exposure studies are widely used as a model to assess the interactions between air pollutants and cells
- They are more economical than animal studies
- Responses either in the cells or in the supernatant provide fast information regarding the toxicological properties of ambient aerosols
- Conventional in vitro assays cannot assess accurately the toxicity of aerosols as exposure studies are limited by the sample collection system:
 1. loss of physical and chemical integrity of the particles during collection
 2. transformations during extraction and resuspension

OUR APPROACH

To eliminate the need for filter-base collection, and reduce sampling artifacts we have developed two new devices that collect ambient aerosols directly:

- 1) into liquid as concentrated suspension
- 2) into cell culture systems

These systems use the three-stage, moderated laminar flow water condensation technology developed by Hering et al. (2014) to grow and collect particles down to 6 nm in liquid



AEROSOL INTO SUSPENSION COLLECTORS



Advantages of the ASCs:

1. Direct collection into liquid reduces sampling artifacts
2. No need for extraction and resuspension steps
3. Collects soluble particle-phase components, and insoluble particles suspended in the collection liquid
4. The concentrated nature of the collection allows reducing sampling times
5. Provides concentrated aerosol/particle samples ready to be characterized using in-vitro and in-vivo assays

1. AEROSOL INTO SUSPENSION COLLECTOR (ASC)

Specifications:

1.5, 3 or 5 lpm sampling flow rate
200-500 μL collection volume



Collection vial



2. DIRECT *IN-VITRO* EXPOSURE SYSTEM (DIVE)


Specifications:

8 lpm sampling flow rate

32 impaction nozzles to minimize disruption of the cell culture medium

CO₂ introduced to the incoming flow (5%)


Temperature (37°C) and RH (>90%) maintained in the collection chamber



Impaction Nozzles




Sample delivery



COLLECTION EFFICIENCY

Collection efficiencies **>90%** for particles sizes between **5 nm** and **10 μm** at moderated temperatures




FIELD DEPLOYMENT: Michigan State University

a) Indoor (Lab testing) - March 2014

- Biosafety level 2 (BSL-2) lab - HEPA filter
- Low #/cc ($\sim 4 \times 10^3$) – lab ambient ai

b) Ambient - May 6-9, 2014

- Empty office room on the first floor
- 100-150m back side of the parking lot
- 3 hr samples, morning and afternoon



FIELD DEPLOYMENT: *In-vitro* assays

Cells:

- Mouse leukaemic monocyte macrophage cell line (**RAW 264.7**)
- Human Bronchial Epithelial cell line (**BEAS-2B**)
- Normal Human Bronchial /Tracheal Epithelial cell line (**NHBE**)

Cytotoxicity:

- Trypan Blue
- MTS Assay


Stimulation:

- **Plating:** 30,000 cells/well, 24-well plates, 1 mL BEGM, 24 hr
- **Stimulation:** **0.001 m³** air/mL (incense burning), 300 μ L/well, 16hr
0.05 m³ air/mL (outdoor), 300 μ L/well, 16 hr

Endpoints:


- **Pro-inflammatory markers:** interleukin 6 (IL-6)
interleukin 8 (IL-8)

RESULTS: *Cell viability*



ASC


Exposure to 0.05 m³/mL of outdoor PM did not have a significant impact on BEAS-2B cell viability



DIVE

No major differences in viability were observed between the control and exposed cells

RESULTS: *Pro-inflammatory effects of PM (ASC)*




➤ Ambient PM induced a:
10- to 34-fold increase in IL-6
5- to 25-fold increase in IL-8


➤ Important diurnal and daily variability in the pro-inflammatory capacity of PM

➤ Higher effect observed in morning periods

RESULTS: *Pro-inflammatory effects of PM (DIVE)*



- Differences between control and exposed cells were observed
- Higher difference for ambient PM



- A significant increase in IL-8 production was only observed in 2 of the days
- Different behavior for different days

SUMMARY

- A. We have developed a new collection system for collecting ambient aerosols directly into highly concentrated liquid suspensions (ASC)
- B. DIVE prototype provides the means to assess aerosol toxicity using *in-vitro* assays under realistic ambient and physiological conditions
- C. Collection Efficiencies for both systems were higher than 90% for particles sizes ranging from 5 nm to 10 μm
- D. Cell viability after exposure to ambient PM was similar to control
- E. Inflammatory responses measured as IL-6 and IL-8 production were observed for exposures to ambient particles down to 3 hrs

ACKNOWLEDGEMENTS

Paula Palmiter for providing the space for ambient collection

Dr. Ryan Lewandowski for his help during field deployment

Grants: 1 RC3 ES019081-01 & 1R43-ES022523

For Product information, contact **Pat Keady** at:

Aerosol Devices Inc.

info@aerosoldevices.com

Website: aerosoldevices.com

ありがとう